

Preparation of Probe for *In Situ* Hybridization

(Rivera lab)

1. Prepare reaction mix

50 ng/ μ l	1.0 μ l	1 μ g/ μ l DNA template (Linear plasmid, PCR fragment)
1X	2.0 μ l	10X Transcription buffer
1X	2.0 μ l	10X DIG RNA labeling mix (Roche Cat # 1277073)
0.01 M	2.5 μ l	0.1M DTT
1 U/ μ l	0.5 μ l	40 U/ μ l RNase inhibitor (RNAsin, Promega N2611)
2 U/ μ l	2.0 μ l	T3, T7 or SP6 RNA Polymerase (20 U/ μ l)
	10.0 μ l	depc ddH ₂ O

2. Incubate at 37 °C for 2 hours.
3. Stop the reaction by adding 1 μ l of 0.5 M EDTA.
4. Add 2.5 μ l of 4 M LiCl and 75 μ l of cold 100% ethanol, to precipitate.
5. Chill at -20 °C for 2 hours (or 30 minutes at -60 °C or below)
6. Centrifuge at 13,000g for 5 minutes.
7. Wash the pellet with 70% ethanol and let dry.
8. Resuspend in 200 μ l of depc TE and add 1 μ l of RNase inhibitor.
9. Check the RNA probe by running 1 μ l on an agarose gel. The signal from the RNA should be 10X stronger than that of the DNA template.
10. Prepare aliquots of 20 μ l and store at -20 °C. It can last for 1 year.
Use ~ 20 μ l (0.1 - 1 μ g) per wholemount in situ hybridization assay.